



Amari Giacomo 10698309

Milazzo Alessandro Giuseppe 10766415

Traversi Gloria 10816956

Usai Chiara 11005835

Zanoni Chiara 10784594



Focus group Proposal Conclusion



'ci ha restituito un pezzo di vita'



ANDREA

Experience: various races

Issues: safety, comfort, upper limbs, outdoor

Purpose: use with families and friends



DAVIDE



ALESSANDRO



ANDREA

Experience: various races

Issues: safety, comfort, upper limbs, outdoor

Purpose: use with families and friends



DAVIDE

Experience: various races

Issues: transporting the trike, customization, comfort, safety, outdoor

Purpose: autonomous training



ALESSANDRO



ANDREA

Experience: various races

Issues: safety, comfort, upper limbs, outdoor

Purpose: use with families and friends



DAVIDE

Experience: various races

Issues: transporting the trike, customization, comfort, safety, outdoor

Purpose: autonomous training



ALESSANDRO

Experience: training

Issues: safety, comfort, upper limbs, outdoor

Purpose: total body training

Purpose: use with families and friends

Purpose: autonomous training

Purpose: total body training

Critical issues





FUNCTIONAL REQUIREMENT

Description of the requirement:

- 1. Prevent further pressure sores
- 2. Be comfortable
- 3. Prevent unwanted movements (of the rider and the device)
- 4. Be lightweight and NOT bulky
- 5. Adaptable to the trike

The independent cell cushion is designed to adjust the air within the chambers in response to pressure changes detected by the MAT, thereby reducing overload and enhancing load balancing.



CUSHION



PRESSURE SENSITIVE – MAT





CUSHION

Adjustable cushion with separate air chambers for personalized support



PRESSURE SENSITIVE – MAT





CUSHION


Material selection: hypoallergenic, breathable, non-slip, waterproof

Design optimization



PRESSURE SENSITIVE – MAT








PRESSURE SENSITIVE – MAT

Sensor-equipped mat that monitors pressure distribution in real time









PRESSURE SENSITIVE - MAT

Sensor & system configuration: type, number, positioning, techincal specifications











MICROCONTROLLER

Integrated circuit that automatically regulates the pressure in the air chambers







MICROCONTROLLER

System configuration & control: hardware and software

